

VIVE PRECISION SPHYGMOMANOMETER

DMD1013BLK

HCPCS: A4663

Track systolic blood pressure measurements conveniently and accurately at home or on the go with Vive Precision's aneroid sphygmomanometer which includes a convenient zipper pouch for storage. The large dial is easy to read and the included D-ring allows individuals to easily secure the cuff for accurate self measurements.

SPECIFICATIONS

What's Included:

- ♥ One aneroid sphygmomanometer
- ♥ D-ring
- ♥ Zippered travel case
- ♥ Manual

Color(s): Gray cuff with black piping

Care:

- ♥ All parts can be washed with a cotton cloth and water.
- ♥ Hand wash nylon cuff in warm water and lay flat to dry.
- ♥ Do not fully submerge the manometer unit.





RELIABLE AND ACCURATE READINGS

Confidently achieve reliable readings and track blood pressure levels with an accuracy of ± 3 mmHg. The calibrated adult aneroid sphygmomanometer measures systolic blood pressure ranging from 0-300 mmHg. Easily portable, a large zippered travel bag provides safe storage at home or on the go.



CONVENIENT SELF MEASUREMENTS

Secured with durable fastener, the lightweight nylon cuff also includes a D-ring to ensure a proper fit for most individuals. Additionally, the included inflation bulb and spring air-release valve enables easy and accurate self-measurements.

EASY TO READ DIAL

Aluminum alloy manometer is stick and crack resistant. The manometer uses a large, high-contrast dial that is easy to read. Perfect for seniors and those with low vision.



LATEX-FREE FLEXIBLE CUFF

Adult sized nylon cuff fits arm circumferences ranging from 9.84" to 15.75". The cuff is easily secured in place with premium velcro to ensure accurate results. The double tube PVC inflation bag has a large spring air-release valve and tapered end valve.

CALL NOW!



www.vivenow.health



(239) 228-8939

CALL NOW!

SR24 - 9R1



VIVE

INSULIN TRAVEL CASE

DMD1040

Providing safe, cool storage for insulin and other perishable medications, the protective Vive cooling case is lightweight and compact. The insulated hard shell case includes two reusable ice packs, a removable organizer and a convenient carry handle. Designed to store insulin needles, pens, refills and other supplies, the organizer features three pockets, a pen holder and internal thermometer. Easily fitting in a purse, backpack, tote or carry on, the cooling case is TSA compliant and insulated to maintain a cool temperature for up to eight hours.

SPECIFICATIONS

Dimensions:

➤ 8.5" by 4.5"

What's Included:

- Vive Precision Cooling Case
- Two ice packs

Care Instructions:

- Remove ice packs and rinse under cool water.
- Dry immediately with a lint-free cloth.
- Spot clean the case and folio with a damp cloth and mild detergent.
- Dry immediately with a lint-free cloth.

Materials:

- Hardshell EVA, Polyester, Latex-free ice packs.

Product Size:

- One size

Call now!

or send us an e-mail
now to find out more!



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com



INSULATED COOLING STORAGE

Safely travel with perishable medications for up to eight hours with the Vive cooling case. Featuring multiple layers for effective insulation, the hard shell travel case also utilizes a tight sealed interior design to aid in maintaining a proper temperature.



INCLUDES TWO ICE PACKS

Perfectly sized for the protective cooling case, the set includes two, water-filled ice packs. The latex-free ice packs are reusable, retaining a cooling temperature for up to eight hours.



REMOVABLE ORGANIZER POUCH

Holding a variety of small items, the cooling case includes a removable folio to safely and neatly store necessary supplies. The pouch includes an interior Fahrenheit thermometer, three pockets and one pen holder to securely hold up to two insulin needles, three pens and two pen refills. The organizer pouch is secured to the case with two strong fastener strips.

PROTECTIVE TRAVEL CASE

Constructed with a durable EVA material, the hard shell case protects fragile contents during regular use. Lightweight and compact for easy portability, the protective cooler case contains two mesh pockets for securing the included ice packs in place. The compact case also features a small carry handle and dual zippers for easily accessing the contents.



(239) 220-5772



sales@vivehealth.com

vive
health

SR25- 3R1

VIVE PRECISION PULSE OXIMETER

DMD1066

HCPCS: E0445

Monitor your blood oxygen saturation with the fast, accurate, and out-of-the-box ready Vive Pulse Oximeter. With clinical accuracy, you can measure your SpO2 and pulse and better track your health. The lightweight design is easy to use with one-touch operation, and easy to read with the bright OLED display.



SPECIFICATIONS

What's Included:

- ♥ Pulse Oximeter (1x)
- ♥ AAA Batteries (2x)
- ♥ Lanyard (1x)
- ♥ Case (1x)

Color(s):

- ♥ Black - DMD1066BLK
- ♥ White - DMD1066WHT (Case and Lanyard not included)

Size(s): One size

Product Dimensions: 2.36"x1.46"x1.26"

Latex free: Yes





READY TO USE OUT OF THE BOX – COMPLETE KIT INCLUDED

The Vive Pulse Oximeter comes fully equipped and ready for use, featuring two AAA batteries and a convenient lanyard for easy portability. Whether you're at home or on the go, you can start monitoring your health right away. The included accessories make it simple to carry and use hands-free, perfect for travel or daily use.



RELIABLE CLINICAL PRECISION AND SPEED

Rely on the SpO2 and pulse readings of this pulse oximeter, as it guarantees accuracy and offers swift results. Its user-friendly design allows for one-touch operation, ensuring efficiency and convenience.

CLEAR RESULTS ON OLED DISPLAY

Experience strain-free result readings on the bright OLED display. The large numbers provide easy visibility, presenting SpO2 and pulse readings clearly and efficiently.



EFFORTLESS MOBILITY WITH LIGHTWEIGHT, COMPACT DESIGN

Ideal for on-the-go health monitoring, this pulse oximeter is lightweight and compact, making it travel-friendly. Safely store it in the included travel case when not in use. Its slim design fits comfortably in a purse or travel bag, ensuring hassle-free portability.

CALL NOW!



www.vivenow.health



(239) 228-8939

CALL NOW!

SR24 - 11R5

VIVE

SILLA PARA DUCHA

LVA1009
HCPCS: E0240

Un elemento esencial para la ducha para aquellos que tienen dificultades para estar de pie durante períodos prolongados o problemas de equilibrio. El asiento extra ancho está equipado con asas para facilitar el proceso de levantarse y sentarse lo más fácil posible. La base se ajusta a la mayoría de las bañeras y regaderas que miden 23.3" de ancho y está equipada con tapas con punta de goma para evitar resbalones durante su uso.



ESPECIFICACIONES

Materiales:

- ▷ Aluminio
- ▷ Material compuesto
- ▷ Espuma de celda cerrada.

Instrucciones de cuidado:

- ▷ Limpiar con un paño húmedo y jabón suave o desinfectante.
- ▷ Secar al aire.

Dimensiones del producto:

- ▷ Asiento: 15.75 x 13"
- ▷ Base: 17.75 x 17"
- ▷ Altura ajustable del asiento: de 16 a 21"
- ▷ Altura del reposabrazos: 6.5" por encima del asiento.
- ▷ Altura del espaldar del asiento: 12.5" por encima del asiento.

Incluye:

- ▷ Una silla para ducha
- ▷ Garantía de por vida.

Call now!

or send us an e-mail
now to find out more!



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com



BASE COMPACTA ANTIDESLIZANTE

Equipada con tapones de goma antideslizantes para mayor estabilidad y seguridad.



ENSAMBLAJE FÁCIL Y SIN HERRAMIENTAS

Configuración sencilla para una personalización fácil.



DISEÑO CONTORNEADO

Múltiples agujeros de drenaje en el asiento eliminan el riesgo de deslizamiento.

DUCHA SEGURA

Siéntete cómodo y seguro al ducharte. Diseñado ergonómicamente, el asiento extra ancho soporta hasta 300 libras. Los mangos acolchados en cada lado proporcionan un agarre cómodo sin riesgo de resbalar. Además, cada pata de la silla está cubierta con una punta de goma antideslizante para mayor seguridad.

DURABLE Y VERSÁTIL

La construcción de aluminio resistente a la corrosión es la elección perfecta para una silla de baño ligera o un banco de ducha que no se oxidará. Ajusta fácilmente la altura del asiento de la silla de ducha para crear una experiencia segura y personalizada. Con un rango de 14 a 19.5 pulgadas, el asiento se ajusta para adaptarse a cualquier altura o para proporcionar estabilidad adicional en pisos irregulares. El espaldar también es removible si se desea.

ASIENTO ANTIDESLIZANTE CÓMODO CON DRENAJE

Diseñado para ofrecer una comodidad superior, el asiento y el respaldo antideslizantes están fabricados con un material compuesto ligeramente texturizado. Los orificios de drenaje reducen el deslizamiento en el agua acumulada, creando un asiento de ducha seguro y cómodo.

MONTAJE FÁCIL SIN NECESIDAD DE HERRAMIENTAS

Ensambla fácilmente el asiento de ducha ligero sin necesidad de utilizar herramientas. La extracción sin herramientas de las patas, reposabrazos y espaldar facilita el almacenamiento y el transporte de manera sencilla y conveniente.



(239) 220-5772



sales@vivehealth.com

**vive
health**

SR24 - 3RI

TRANSFER BENCH

LVA1041

vive
vivehealth.com




To see all FAQs
in one place,
visit vhealth.link/30d1d

OWNER'S MANUAL

TRANSFER BENCH

WHAT'S INCLUDED



- 1 x Small Seat Section (A)
- 1 x Large Seat Section (B)
- 2 x Bench Legs w/ preinstalled Suction Cup Feet (C)
- 2 x Bench Legs w/ preinstalled Capped Feet (D)
- 1 x Backrest Assembly (E)
- 1 x Arm Rest (F)




ASSEMBLY INSTRUCTIONS



To see a video demonstration, check out vhealth.link/041

1. Turn the Seat Sections (A) and (B) over so the leg sockets are facing up.
2. Connect the two male frame ends on the Small Seat Section (A) into the female frame ends under the Large Seat Section (B) by depressing the locking pins and pushing the seat sections together until the pins pop into place with an audible click.
3. Fit each Bench Leg w/ preinstalled Suction Cup Foot (C) into the leg sockets on the Small Seat Section (A). Make sure each leg is fitted in place, then align the holes in the Leg with the holes in the socket and secure the pieces in place with the locking pin.



NOTE - The plastic stoppers fitted at the end of each Leg do not need to be removed before installation.

- Repeat Step 2 with the other Bench Legs (D) in the sockets on the Large Seat Section (B).

If needed, adjust the height of each leg by depressing the buttons at the bottom of each leg. With the button depressed, slide the leg up or down until the pin aligns with an adjustment hole at the desired height. The pin will click into place. Make sure the legs are adjusted and the seat is level when installed in the bathtub.

- Turn the bench onto its long side and attach the Backrest Assembly by depressing the locking buttons and fitting the backrest bars into the tubes on the underside of the bench. The backrest can be attached on either side, to make the bench reversible depending on the configuration of your shower.



- Stand the bench up on its 4 feet and attach the arm rest by depressing the locking buttons and fitting the bar into the sockets on one side of the bench. Make sure the armrest bar curves away from the seat.



WARNINGS

- Seat back is not height adjustable.
- 400 lbs weight limit.
- Wipe down with dry cloth after use.
- Make sure all components are securely connected before sitting or leaning on the bench.
- Do not sit on the bench if the seating surface is not level.
- Check periodically to ensure that water is not collecting inside the legs. If water collects inside, remove the legs, drain the water inside, and reattach.

HAVE MORE QUESTIONS?

Check out our list of Frequently Asked Questions at
[vhealth.link/30d1d](https://vivehealth.com/health.link/30d1d) for helpful answers.



And if that doesn't answer your question, our
customer service team would love to help!
Feel free to connect with them by phone,
e-mail, or chat on our website

service@vivehealth.com
1-800-487-3808
vivehealth.com

Distributed by

vive
health

8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com



VIVE FOLDING CANE

MOB1001

HCPCS: E0100

Perfect for active individuals, the Vive Folding Cane is lightweight and compact, quickly collapsing to fit into a purse, tote or the included storage pouch. Supporting up to 250 pounds, the durable folding cane locks into place for stable use on any surface and is capped with a nonslip, rubber cane tip for greater traction on slippery surfaces. Fitting comfortably in either hand, the foldable cane is topped with an ergonomic, molded handle for a secure grip that minimizes hand and wrist fatigue. Height adjustable, the cane stands 33" - 37" tall and includes a locking storage clip and travel pouch for added convenience.

SPECIFICATIONS

Materials:

➤ Anodized aluminum, Rubber tip and grip

What's included:

➤ Folding Cane

Dimensions:

➤ Height: Adjustable at 1" increments ranging from 33"-37"

➤ Diameter: 3/4"

Color:

- Black
- Blue
- Teal
- Purple
- Dark Butterfly
- Leopard
- Blue Burst
- Grey Poly

Call now!

or send us an e-mail
now to find out more!



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com



FOLDING DESIGN

Convenient for storage or travel.



EASILY PORTABLE

Slips into your purse, briefcase, travel bag, wheelchair bag and more.



ERGONOMIC GRIP

Contours comfortably to your hand to reduce fatigue.

COLLAPSIBLE AND COMPACT

Providing a convenient option for those needing assistance, the Vive Folding Cane smoothly and securely locks in place to provide extra stability while walking on any surface. Folding at four points along the cane shaft, the collapsible cane is compact, easily fitting into a purse, briefcase, carry on or wheelchair bag.

LIGHTWEIGHT ADJUSTABLE ALUMINUM FRAME

Constructed with a durable, corrosion-resistant anodized aluminum, the Vive Folding Cane is exceptionally lightweight yet sturdy, safely supporting up to 250 pounds. The cane is easily adjusted in one-inch increments from a minimum of 33" to a maximum height of 37".

ERGONOMIC HAND GRIP AND STRAP

Ergonomic hand grip contours to the hand to reduce hand fatigue, prevent muscle cramping and relieve pressure on the wrist. The integrated wrist strap keeps the cane from falling to the ground when released.

SLIP-RESISTANT RUBBER TIP

Collapsible cane is capped with a slip-resistant rubber tip for additional traction and stability. The cane is also compatible with most replacement tips, including tripod tips, quad tips and ice tips.



(239) 220-5772



sales@vivehealth.com

vive
health

SR24 - SR1

vive health

portal.vivehealth.com

VIVE ROLLATOR WALKER

MOB1010BLK

HCPCS: E0143

Oversized wheels with a dual braking system ensure stability on any surface, while easily accessible storage conveniently holds wallets, phones, drinks and more. The compact rollator folds down to easily fit in tight spaces with its hinged and removable wheels.

WHAT'S INCLUDED

- Rollator Frame
- 2 Removable Front Wheels
- 2 Fixed Back Wheels
- Removable Storage Bag

COLORS 



SAFELY AND SECURELY MANEUVER INDOORS AND OUT



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com

Specifications

Overall Dimensions 25" x 28"

Folded Dimensions 10" x 20" x 32"

Seat 17.5" x 9.5"

Wheels Diameter 8"

Storage Bag Dimensions 15" x 13" x 3.5"

Storage Bag Weight Capacity 11 lbs

Adjustable Handle Height 31" to 36"

Handle Height can be adjusted in 1" increments with the push pins located on the sides of the rollator

Materials Lightweight Aluminum

Product Weight 19 lbs.

Weight Capacity 300 lbs.

SUPERIOR SUPPORT AND STABILITY

Lightweight and durable, the Vive rollator provides stable support, and easily maneuvers over any type of surface with large 8" sport wheels. The handle height is adjustable for a customized fit for every individual. Supporting up to 300 pounds, the lightweight aluminum frame is corrosion-resistant and folds down for easy storage and transport. Includes a wide seat and backrest for comfortably resting while on the go.

COMPACT FOLDABLE DESIGN

The collapsible rollator is compact when folded, easily fitting in the trunk of a car or even in the backseat, allowing for hassle-free folding in one easy step. Simply pull up on the center strap under the seat. For more compact storage, the front wheels are removable with a simple push button release and the back wheels are hinged to fold up against the frame.

DUAL BRAKING SYSTEM

Each of the ergonomic handles is angled to position hands in a neutral position, reducing fatigue and allowing easy access to the bicycle-style hand brakes. Wheel locking brakes ensure the rollator will not move when activated.

CONVENIENT SAFE STORAGE OPTIONS

Conveniently and safely store daily necessities in the detachable rollator bag. Safely holds a purse or wallet, medications, electronic devices and more. Also includes an integrated cane holder.





VIVE

HEAVY DUTY WHEELCHAIR

MOB1078

Experience a durable yet lightweight mobility solution designed for a variety of individual needs. With a strong steel frame, this wheelchair can support up to 400 pounds. It can also be folded for easy transportation and storage. The 20-inch wide seat and padded armrests ensure comfort, while the adjustable footrests allow for a more personalized fit.


WHAT'S INCLUDED

- 1 Heavy Duty Wheelchair

COLORS:



HEAVY-DUTY AND LIGHTWEIGHT



Product Dimensions:

Seat width:	20"
Seat height:	19.50"
Seat depth:	16"
Back height:	18"
Total width:	27"
Total height:	36.3"
Total length(with footrest):	42"
Max loading:	400LBS



HEAVY-DUTY AND LIGHTWEIGHT

The Vive Foldable Wheelchair Transport offers robust support, accommodating up to 400 pounds. With a lightweight steel frame, this wheelchair is portable and makes for an excellent travel companion. It also folds, for hassle-free storage when not in use.

ADJUSTABLE FOR GREATER COMFORT

Experience superior comfort with this wheelchair's 20-inch wide seat, providing extra space for a more relaxed sitting position. You can customize your comfort further with the adjustable footrests, which conveniently swing away if needed. The padded armrests ensure comfort, while the easy-to-clean material of the chair ensures a cozy and hygienic experience.

ENHANCED MANEUVERABILITY

Navigate easily thanks to the effortless push design and the durable wheels on the Vive Foldable Wheelchair Transport. The hand-breaking system enhances control and safety, for greater peace of mind when maneuvering the chair. This wheelchair fits through most doorways with a total width of 27 inches, so you can move between indoor and outdoor spaces.

VERSATILE DESIGN FOR A VARIETY OF NEEDS

With four sturdy wheels and an adjustable footrest, the Vive Foldable Wheelchair Transport is great for individuals with a variety of needs. Provide seamless transportation between rooms and from hospitals to your home, or use this wheelchair to add stability during travel. With footrests that swing open, this wheelchair is easy to access and exit, making it great for those who are recovering from surgery or injury.



VIVE
LIGHTWEIGHT ROLLATOR
MOB1079

Crafted with precision and weighing only 14 pounds, the lightweight aluminum frame ensures effortless transportation and storage when not in use while safely accommodating up to 300 pounds. The ergonomic dual-braking system provides enhanced safety with secure wheel locking when not in use. Enjoy superior comfort and support with the padded seat and backrest, complete with a handy storage compartment for your essentials. Whether you are navigating indoor spaces or exploring the outdoors, the 6-inch wheels guarantee smooth movement, making it an ideal companion for travel or everyday use.

WHAT'S INCLUDED

- 1 Lightweight Rollator

COLORS ■ ■

14 lbs
LIGHTWEIGHT





portal.vivehealth.com

Specifications

Width	22" x 24"
Adjustable Height	31" to 36"
Seat	14 x 12"
Seat height from the floor	21.5"
Wheels	6"
Materials	Aluminum
Product Weight:	14lbs
Weight Capacity	up to 300 pounds

LIGHTWEIGHT ALUMINUM FRAME

The carefully crafted, lightweight aluminum frame of the Lightweight Rollator weighs a mere 14 pounds, allowing you to easily transport and store this rollator when it is not in use. The strength of the aluminum frame safely supports up to 300 pounds, for greater peace of mind.

ERGONOMIC DUAL-BRAKING SYSTEM

With comfortable handles that allow a natural grip and relieve tension on your wrists and hands, this rollator provides an ergonomic braking system that is easy to use. With dual braking, the wheels lock securely in place when not in use, for enhanced safety.

SUPERIOR COMFORT AND SUPPORT

Sit with greater comfort thanks to the padded seat and backrest. The integrated storage compartment allows you to keep your essentials close at hand, for greater convenience while on the go.

SMOOTH MOVEMENT INDOORS AND OUTDOORS

The 6" wheels are crafted to seamlessly glide over any surface, whether you are exploring new outdoor spaces or moving about a familiar indoor area. Take this rollator walker with you while traveling or use it in the comfort of your own home.

vive health 8955 Fontana Del Sol Way
 Naples, FL 34109
 1-800-487-3808
www.vivehealth.com

SR23-12R1



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com



VIVE

CORE SHOWER COMMODE TRANSPORT WHEELCHAIR MOB1115

The Core Shower Commode Transport Wheelchair is a compact and versatile 3-in-1 mobility solution, effortlessly functioning as a rolling chair, shower seat, and commode. Designed with a lightweight, rust-resistant aluminum frame, it ensures durability without unnecessary bulk. Locking wheels provide stability, while the comfortable seat, supportive backrest, and foldable footrest enhance user comfort. Its space-saving, foldable design allows for easy storage and transport, making it ideal for those seeking convenience and reliability.

SPECIFICATIONS

What's Included:

- Rolling chair and bucket

Product Dimensions:

- (Lowest and Highest possible): min: 36" max: 38"
- Base Dims (Unfolded): 22" x 28.5"
- Front Wheel Size: 4"
- Back Wheel Size: 4"
- Dims when folded: 22" x 20.5"
- Seat Height (if there is a seat): 21"
- Seat Dims: 17" x 17"
- Backrest Size: 18" x 6.75"
- Width between handles (Applies if there is a seat): 17.5"
- Bucket Dims: 12" x 9" x 6"
- Foldable Footrest: 21" x 10"

Product Material(s): Aluminum and plastic

Call now!

or send us an e-mail
now to find out more!



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com



VERSATILE 3-IN-1 DESIGN

Easily transitions between a transport chair, shower seat, and commode, making it a practical solution for individuals with mobility challenges.



LIGHTWEIGHT YET DURABLE

Constructed from rust-resistant aluminum, this wheelchair offers exceptional durability while remaining easy to maneuver and transport.



ENHANCED STABILITY & COMFORT

Locking wheels ensure secure positioning, while the padded seat, supportive backrest, and foldable footrest provide added comfort and ease of use.

COMPACT & TRAVEL-FRIENDLY

Designed for convenience, the foldable frame allows for effortless storage and portability, making it ideal for home use or travel.

COMPLETE & READY TO USE

Comes with a rolling chair and bucket, providing everything needed for a comfortable and hassle-free experience.



(239) 220-5772



sales@vivehealth.com

vive health

SR25-2R1

CORE TRANSPORT WHEELCHAIR

MOB1131

vive®
MOBILITY



To see all FAQs in
one place, visit
vivehealth.link/be6ca1

vivehealth.com


CORE TRANSPORT WHEELCHAIR

FOLDING / UNFOLDING THE CHAIR



For a video demonstration, check out vhealth.link/3074f2



1. From a stable position, grip the armrests or the chair frame under the sides of the seat and pull the frame away from the center of the chair until the fabric seat is fully unfolded.



2. Grab the folded down hand grips and pull them up until the backrest upholstery is full extended and both hand grip bars snap into place locking them at a 90 degree angle.



3. To fold the chair, Push the two black tabs beneath the hand grips up and in towards the frame to disengage the locking pins. You can then lower the handles down. Next, lift the chair by the seat strap until the seat folds together, or lift by the arm rests and push the two sides together until the seat fold together.



ATTACHING / REMOVING THE FOOT RESTS



For a video demonstration, check out vhealth.link/4b8926

1. Fit the foot rest over and down onto the two metal placement pins as shown.



2. Then rotate each foot rest forward until the silver clasp locks into place and the foot rest is securely fastened to the chair frame.



3. To remove the foot rests, pull the wide end of the silver clasp in towards the middle of the chair. You will feel the clasp release its grip on the frame, and the foot rest will swing out where it can be lifted off the two metal placement pins.

ADJUSTING THE HEIGHT OF THE FOOT RESTS

Pull up on the black locking tab on the back to unlock, then adjust the footrest height by pressing in the locking pin on the pipe of the footrests. Adjust the height up or down as needed making sure the pin depresses back into the hole for the height you need. Lock the black tab on the back when satisfied with your height.



LOCKING THE HANDBRAKES

The parking brake is used to help prevent the Wheelchair from moving when not in use, or to prevent it from moving while sitting or standing. To engage the parking brake, push the parking brake handle on each on the rear legs forwards toward the front of the wheel. To release the brake pull the same handle backward on each side.



MAINTENANCE

Tire Maintenance

Tires should be cleaned occasionally with a damp cloth. Replace if tires become severely worn or cracking appears. All casters should be checked periodically for cracking and lubricated occasionally.

General Maintenance

Wipe metal parts with a soft cloth at least once a week. All aluminum parts can be cleaned with auto wax or similar products. If exposed to moisture, dry immediately. Do not use an abrasive cleanser that will scratch the finish.

CAUTION

Please note that this chair is designed to use for transporting another individual and is not meant for self transportation. This Transport Chair is designed to safely be a stable center of gravity but it is important to try and keep a stable center of gravity of all times while performing activities that may change your center of gravity or a passengers center of gravity. Never reach out so far that is would require you to shift in the seat or bend over out of the chair.

When you have to reach for objects that require you to lean, be certain your front wheels are pointing forward as this will stabilize the chair.

Ramps

Be sure of your own capabilities, limitations and strength before attempting to push another individual up a ramp. Correct techniques on ramp negotiation should be shown to you by a qualified professional. If you must stop on an incline, avoid any and all sudden or abrupt movement when you resume movement.

Because of backward angle of the chair, sudden abrupt movement could cause instability. Do not attempt to change direction when moving in incline or decline, this could result in tipping.

WARNING 

- Wheelchair will support up to a 250 lb individual. Do not exceed chair weight limit.
- Do not attempt to use wheelchair without reading carefully all instructions in this manual.
- This wheelchair has not been approved as a seat that can be used within vehicles of any kind.
- Always transfer wheelchair user to a vehicle manufacturers approved seat and use safety restraints available.
- Do not operate this wheelchair on roadways, streets or any surfaces with vehicular traffic.
- Do not attempt to operate chair on inclines without proper assistance.
- Wheelchair should be on a smooth, stable and level area with the hand brakes engaged by the transporter before transferring the passenger in or out of the wheelchair.
- When transferring in or out of wheelchair, do not stand on or apply weight to foot rests.
- Do not lean over or apply weight to the top of wheelchair back. This can cause wheelchair to tip over.
- Do not stand on wheelchair.
- Abruptly changing directions while going down an incline could cause instability and/or increase the possibility of tipping over.

HAVE MORE QUESTIONS?

Check out our list of Frequently Asked Questions at
[vhealth.link/be6ca1](https://vivehealth.com/health.link/be6ca1) for helpful answers.



And if that doesn't answer your question, our
customer service team would love to help!
Feel free to connect with them by phone,
e-mail, or chat on our website

service@vivehealth.com
1-800-487-3808
vivehealth.com

Distributed by

vive
health

8955 Fontana Del Sol Way
Naples, FL 34109
Made in China

FOLDING PEDAL EXERCISER

QUICK START GUIDE

OVERVIEW

We are constantly answering questions and recording helpful videos to make using your Vive Folding Pedal Exerciser as easy as possible. Check out the included link and QR code to help you through the process.



To see all FAQs in one place, visit vhealth.link/r14

WHAT'S IN THE PACKAGE

- 1x Folding Pedal Exerciser

QUICK SETUP

1. Unfold the Pedal Exerciser legs. You will hear an audible click when the legs securely lock into place.
2. Press and hold the red Mode button to turn on the screen.

3. Place the Pedal Exerciser on a flat surface.
4. Begin pedaling, using either your arms or legs. The screen will begin counting.
5. Turn the adjustment knob to increase or decrease the pedal tension.
6. When you're done working out, the screen will automatically shut off after a few minutes of inactivity.

HOW TO USE THE DISPLAY SCREEN



Component	Description
Red (Mode) Button	Turns the Folding Pedal Exerciser on. Hold the button down to reset any previous information to zero. Press it quickly to cycle through the four (4) different display modes: Time, Calories, Revolutions per Minute, and Scan.

CNT (Count) Display	Shows the current number of full revolutions on the Pedal Exerciser since the last reset.
TIME Indicator	Shows when the Exerciser is in Time Display mode, showing how long the current workout has lasted. This is the default display mode.
CAL (Calories) Indicator	Shows when the Exerciser is in Calories Display mode, showing approximately how many calories have been burned during the current workout.
RPM (Revolutions per Minute) Indicator	Shows when the Exerciser is in RPM Display mode, showing approximately how many revolutions are made on average during the current workout.
SCAN Indicator	Shows when the Exerciser is in Scan Display mode, cycling through Time, Calories, and RPM Display modes for the current workout.

REPLACING THE BATTERY

If the display will not turn on, it may require a new battery.

1. Use a flathead screwdriver in the divot on the right side of the display to remove the display from the Exerciser (this may require a bit of force). The display will be wired to the device even after it's popped free, so do not pull forcefully on the removed display.
2. Slide the old battery from the socket on the underside of the display, and replace it with a new LR44 button-cell battery.
3. Insert the display into the left side of the frame in the Exerciser, and press the right side down to snap it into place (this may require a bit of force).



ADDITIONAL INFORMATION AND SAFETY WARNINGS

- If the Pedal Exerciser is wobbling, make sure that it is on a flat, level surface.
- The pedal arms can heat up after extended use. This is normal and is created by the friction in the resistance mechanism. Allow the Exerciser time to rest and cool down before handling to prevent injury.

SPECIFICATIONS

- Height: 9 in.
- Pedal Width: 4 in.
- Depth: 19 in.
- Weight: 5.15 lb.
- Width: 14.5 in.
- Display Battery: LR44 button-cell

vive®

Distributed by

vive
health

8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com
Made in China

VIVE 50-MODE TENS AND EMS

RHB3004

Experience versatile and personalized pain relief with the Vive 50-Mode TENS Unit. This advanced unit offers an impressive array of 50 modes designed to cater to various pain and tension needs. With 20 TENS modes for muscle relaxation and pain relief, 20 EMS modes for muscle strengthening and recovery, and 10 FITNESS modes for targeted muscle group stimulation, you have a comprehensive pain relief solution at your fingertips. The large LCD display ensures effortless navigation through modes, while the Type C charging capability guarantees quick charging and prolonged usage. The package includes a variety of electrode pads, including square, circular, and rectangular shapes to ensure optimal contact for effective relief and relaxation. This comprehensive TENS unit is your go-to for personalized pain management and relief.

SPECIFICATIONS

Materials:

- ✔ Lens:PMMA.
- ✔ PCB:FR4.
- ✔ Buttons, plastic shells: ABS.
- ✔ Electrodes: Glycerine, Water, Black carbon conductive layer, Release paper.
- ✔ Latex-free.

What's included:

- ✔ Tens/Ems Device
- ✔ 1 Type C Charger
- ✔ 4pcs 2 x 2" Square Pads
- ✔ 2pcs 1.4" Circle Pads
- ✔ 2pcs 2x4" Rectangular Pads
- ✔ Plastic Carrying Case

Care Instructions:

- ✔ Device - Wipe clean with damp cloth, Wipe dry.
- ✔ Electrode pads - Wet with water and wipe clean with cloth.

Color:

- ✔ White



Call now!

or send us an e-mail
now to find out more!



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com



EFFECTIVE PAIN RELIEF

Offers an array of pain relief options that can be tailored to help relieve pain and tension.



SPACIOUS LCD DISPLAY

The clear screen makes it effortless to navigate through different modes and adjust settings.



LONG-LASTING BATTERY

The long-lasting battery ensures that you can rely on your device for extended periods.

EFFECTIVE AND PERSONALIZED PAIN RELIEF

The 50-Mode TENS Unit by Vive offers an array of pain relief options that can be tailored to help relieve pain and tension. The first 20 modes, specifically designed for TENS applications, cater to muscle relaxation and pain relief. The 20 EMS modes are tailored for muscle strengthening and recovery, while the 10 FITNESS modes help you achieve your fitness goals by stimulating muscle groups.

LARGE AND EASY-TO-READ LCD DISPLAY

The Vive 50-Mode TENS unit has a spacious LCD display that ensures easy readability. The clear screen makes it effortless to navigate through different modes and adjust settings according to your preference. No need to strain your eyes or struggle to figure out the settings—this TENS unit ensures hassle-free use.

QUICK CHARGE AND LONG-LASTING BATTERY

Charge your TENS unit fast with the hassle-free Type C charging capability. The quick charge feature ensures minimal downtime between sessions, allowing you to focus on managing your pain. Additionally, the long-lasting battery ensures that you can rely on your device for extended periods, whether you are using it at home, at work, or during travel.

INCLUDES A VARIETY OF ELECTRODE PADS

The 50-Mode TENS unit comes complete with a set of 4 square pads, 2 rectangular pads, and 2 circle pads to accommodate different muscle groups and target areas. These electrode pads ensure optimal contact for effective pain relief, muscle stimulation, and relaxation. With these included pads, you have everything you need for a comprehensive pain management experience.



(239) 220-5772



sales@vivehealth.com

vive
health

SR24 - 1R1



VIVE WRIST BRACE

SUP1021

HCPCS: L3908

The Universal Wrist Brace by Vive, ergonomically designed to provide maximum support for the wrist, the Vive wrist brace is great for those with arthritis, carpal tunnel, repetitive stress injury, tendonitis, sprains and strains. The compression brace reduces pain, inflammation and swelling. A removable splint allows you to determine the level of support provided for a customizable experience. Lightweight and breathable, the nonslip brace is suitable for extended wear.

SPECIFICATIONS

Universal Size:

- Fits Left or Right hand

Materials:

- 50% Neoprene
- 30% Nylon
- 15% Polyester
- 5% Spandex

Care and cleaning:

- Remove Splint Before Washing.
- Hand Wash in Cold Water with Mild Soap.
- Hang Dry

What's included:

- 1pc Wrist Brace
- 1pc Splint

Color:



Call now!

or send us an e-mail
now to find out more!



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com



INJURY SUPPORT

Provides support to injured muscles and helps to protect against re-injury



REMOVABLE SPLINT

Allows you to determine the stiffness and support of your brace



EXTRA-STRENGTH FASTENER

Prevents constant readjusting and allows you to customize the fit of the brace

PROTECT AND STABILIZE YOUR WRIST

Stabilizing and supporting the wrist, the lightweight wrist brace by Vive provides maximum support to reduce pain and inflammation and prevent re-injury. Fully adjustable, the compression brace retains therapeutic heat to improve circulation and promote healing. The lightweight brace, with a removable splint, allows the full use of your hand, so you can easily complete everyday activities.

CUSTOMIZABLE COMPRESSION THERAPY

Two adjustable straps allow you to customize the amount of compression to effectively reduce pain, inflammation and swelling. Great support brace for arthritis, carpal tunnel, tendonitis, repetitive stress injury, post-surgery, sprains and strains.

REMOVABLE SPLINT FOR FLEXIBLE USE

A lightweight, durable splint stabilizes the wrist in a neutral position for optimal relief. The splint is also removable, allowing you more flexibility while providing adequate wrist support.

ADJUSTABLE LIGHTWEIGHT DESIGN

Universally sized to accommodate most hand sizes, the wrist support brace features two adjustable straps with extra strength fasteners for a personalized fit without slipping or bunching. The breathable neoprene keeps your hand and wrist cool and comfortable throughout the day.

VIVE GUARANTEE

60 day unconditional guarantee so you can purchase now with confidence.



(239) 220-5772



sales@vivehealth.com

vive
health

SR23-2R1



VIVE HARD NIGHT SPLINT

SUP1035

HCPCS: L4396/L4397

Relieve foot pain and inflammation while you sleep with the Vive Hard Night Splint. Constructed with a hard composite exterior and lined with a soft, breathable padding, the night splint gently stretches the foot, maintaining a neutral position to reduce pain and inflammation related to plantar fasciitis, achilles tendonitis, and drop foot. The adjustable splint includes removable stretch wedges to elevate the foot for maximum pain relief. Use the versatile splint with or without a sock, on either the left or right foot.

SPECIFICATIONS

Sizes:

- Small: Men's: Up to 6.5 - Women's: Up to 7
- Medium: Men's: 7 - 8 - Women's: 7.5 - 9.5
- Large: Men's: 8.5 - 11 - Women's: 10 - 12.5

Materials:

- Breathable Neoprene, Strong fastener & Metal splints.

Care Instructions:

- Hand wash using cool water and mild soap.
- Dry by dabbing with towel to remove the majority of the water.
- Allow insole to air dry until all moisture is removed.

What's included:

- Single hard night splint
- 1 x spike ball
- 60 day guarantee

Call now!

or send us an e-mail
now to find out more!



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com



ALLEVIATE FOOT PAIN AND INFLAMMATION

Splint maintains neutral foot position while creating a gentle stretch.



CUSTOMIZE YOUR STRETCH

Use the wedges separately or combined for your desired stretch.



VERSATILE SPLINT DESIGN

Can be repositioned to fit either foot.

ALLEVIATE FOOT PAIN AND INFLAMMATION

Maintaining a neutral foot position that gently stretches the plantar fascia, the night splint reduces pain and inflammation due to plantar fasciitis, drop foot and Achilles tendonitis. Gently stretching the foot, the splint is worn throughout the night to alleviate painful first steps in the morning. The rigid exterior supports the calf, ankle and foot so you can rest comfortably.

TWO WAYS TO OPTIMIZE THE STRETCH

Customize the amount of stretch for optimal pain relief with the use of two comfortably padded stretch wedges and dual tension straps. The tension straps can be adjusted to increase angle of the foot. Made with a soft foam, the wedges elevate the foot to ensure proper positioning for maximum relief. Use each wedge separately or combined to achieve the perfect amount of stretch.

DURABLE EXTERIOR WITH A SOFT INTERIOR

Durably constructed with a hard composite material, the exterior splint provides exceptional support that can withstand extensive use. The interior of the hard night splint is comfortably padded with a breathable foam with a soft lining. The soft lining prevents skin irritation, blisters and abrasions.

VERSATILE DESIGN

Secured with three wide adjustable straps, the latex-free night splint fits both men and women. The straps are made with a strong webbing and fastener material. Each strap also has a buckle, eliminating the need to adjust the straps every night. Versatilely designed, the splint can be worn on the left or right foot.



(239) 220-5772



sales@vivehealth.com

vive
health

SR23 - BR3



VIVE

ARCH SUPPORTS

SUP1043

HCPCS: L3050

Lifting and supporting the arch, the Vive Arch Support provides adjustable compression to relieve painful foot conditions such as plantar fasciitis, heel spurs, sprains and strains. Including three soft gel inserts, the adjustable compression brace has an integrated pouch, allowing you to choose the right amount of upward compression for your arch. Breathable and lightweight, the low-profile support brace easily fits in most shoe styles and also provides comfortable support for walking barefoot.

SPECIFICATIONS

Sizes:

- Universal Size

Care Instructions:

- Remove gel inserts.
- Machine wash cold.
- Air dry.

Materials:

- 3mm thick neoprene
- Extra-strength fastening material
- Gel support inserts

What's included:

- Universal arch support - one pair
- Six gel inserts
- 60 Day Guarantee

Colors:

- Black
- Grey
- Beige

Call now!

or send us an e-mail
now to find out more!



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com



SUPPORTS PROPER FOOT ALIGNMENT

Lifts upward for greater arch stability.



3 LEVELS OF COMPRESSION

Gel inserts provide the perfect level of support.



BREATHABLE AND LIGHTWEIGHT

Latex-free blend stays cool and dry for exceptional support.

PLANTAR FASCIITIS AND HEEL PAIN RELIEF

Providing relief from painful foot conditions such as plantar fasciitis, heel spurs, strains and sprains, the adjustable compression brace supports the arch while in motion to reduce pain while exercising, hiking or walking. The adjustable arch brace also relieves knee, hip and back pain by maintaining proper foot alignment.

ADJUSTABLE COMPRESSION SUPPORT

Supplying upward compression and support, the arch support brace has an integrated pouch that holds up to three soft gel inserts for added support and comfort. The brace lifts and supports the arch while standing and helps to maintain proper foot alignment during motion, both pronation and supination. Designed for comfort, the low-profile arch brace can be worn with most shoes and sandals, with or without socks or barefoot.

STAYS IN PLACE ALL DAY

Secure the brace in the most comfortable position with the strong fastener strap that also allows the level of compression to be adjusted for optimal support and comfort. The inner lining is made with a non slip material that holds the brace firmly in place, eliminating slipping, bunching and pinching for a comfortable, secure fit.

LIGHTWEIGHT AND BREATHABLE

The lightweight arch support does not add bulk to your foot, easily fitting in most shoes. Constructed with a breathable, latex-free material, the brace also remains cool and dry throughout the day for exceptional comfort.



(239) 220-5772



sales@vivehealth.com

vive
health

SR24 - 4R1



VIVE KNEE BRACE

SUP2009BLK

HCPCS: A4467

Designed with unique spring splint supports, the Vive Knee Brace stabilizes and supports the injured or weak knees during any activity. The flexible knee brace aids in the recovery of sprains and strains, tendonitis and provides relief from arthritis and chronic pain. Reversible for use on the left or right knee, the open patella brace also relieves painful pressure on the kneecap while supporting proper alignment of the ligaments. The latex-free knee brace is constructed with a breathable fabric blend and includes four removable spring splints and two strong fastener straps for a customized fit.

SPECIFICATIONS

- Fit knee circumference: up to 20"
- Above the knee: up to 22.5"
- Below the knee: up to 20.5"

Care Instructions:

- Remove all four spring splints before washing.
- Machine wash in cold water with like colors.
- Air dry completely before replacing the spring splints.

Materials:

- Latex-free nylon
- Spandex blend
- Aluminum splints
- Extra-strength fastening material

What's included:

- Vive Knee Brace
- Four aluminum spring splints

Call now!

or send us an e-mail
now to find out more!



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com



SUPPORTIVE KNEE BRACE

Removable splints provide additional stability.



CUSTOMIZABLE COMPRESSION

Straps easily adjust for desired compression.



OPEN PATELLA DESIGN

Increase range of motion.

STABILIZING KNEE SUPPORT

Providing stabilizing support for weak or injured knees, the Vive Knee Brace relieves pain due to arthritis, sprains, strains, tendonitis and chronic pain. Featuring an open patella, the knee brace supports the knee cap to minimize stress and painful pressure. The knee brace allows for a wide range of motion for daily activities, including exercise, sports and more.

REMOVABLE SPRING SPLINTS

Flexibly designed to support the knee while in motion, the Vive Knee Brace includes four curved aluminum spring supports. Unique stabilizing splints are also gently curved to contour to the knee, providing additional support for proper ACL and MCL alignment. The spring supports are removable for greater flexibility and cleaning.

SLIM LOW-PROFILE DESIGN

Reversible and adjustable, the Vive Knee Brace is great for both men and women, fitting knee circumferences up to 20". The slim, low-profile brace can be worn under regular fit clothing or over clothing as desired.

BREATHABLE FABRIC BLEND

Machine washable for easy cleaning, the supportive knee brace is constructed with a soft, latex-free fabric blend. The knee brace material is perforated for breathability, staying cool and dry throughout the day. Includes two strong fastener straps and a non-slip lining to ensure a comfortable, yet secure fit that will not slide, twist or bunch.



(239) 220-5772



sales@vivehealth.com

vive
health

SR24 - 5R1



VIVE

HAND & WRIST IMMOBILIZER

SUP2024

HCPCS: L3807 / L3809

Stabilizing and supporting the hand and wrist following an injury or stroke, the Vive hand and wrist immobilizer securely holds the palm, fingers, thumb and wrist in a resting position. The immobilizing brace is soft and breathable, effectively relieving pain and inflammation while promoting recovery. Ideal for use both day and night, the stroke splint features an ergonomic rigid splint and four adjustable fastener straps for a secure, yet comfortable fit.

SPECIFICATIONS

➤ Size: Left or Right

➤ XS - length: 9.84"; palm width: 3.74"

➤ S - length: 11.81"; palm width: 3.94"

➤ M - length: 12.60"; palm width: 4.33"

➤ L - length: 13.78"; palm width: 4.53"

Materials:

➤ Composite splint, Lycra/cotton blend, Latex-free foam, fastener material

Care and cleaning:

➤ Remove the rigid splint before cleaning

➤ Machine wash the cover in cold water with like colors

➤ Air dry completely before inserting the rigid splint

What's included:

➤ Vive Hand and Wrist Immobilizer

➤ 60 day guarantee

Call now!

or send us an e-mail
now to find out more!



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com



IMMOBILIZING SUPPORT

Stabilizes the hand and wrist day and night to relieve pain and inflammation due to injury, arthritis and tendonitis



CUSHIONED LINING

Features a soft, velvet cover that is also lightly padded for greater comfort



ADJUSTABLE THUMB LOOP

Separate and support the thumb in a comfortable, neutral position with the adjustable thumb loop

STABILIZING HAND AND WRIST SUPPORT

Ideal for day and nighttime support, the Vive hand and wrist immobilizer holds the affected arm in neutral, resting position to promote healing and pain relief. The brace is softly padded and utilizes a rigid splint to support and stabilize the wrist and hand to relieve pain and inflammation due to injury, tendonitis and arthritis. Also great for supporting proper hand and wrist positioning following a stroke or traumatic injury.

ERGONOMIC SPLINT

Constructed with a rigid composite material, the immobilizing splint is ergonomically shaped for optimal support and comfort. The shaped splint supports a resting hand position to reduce flexion contractures and retain strength.

COMFORTABLY PADDED

Cushioning the arm, wrist and palm, the stroke hand splint is lined with a soft foam padding for greater comfort throughout the day and night. The brace cover is made with a soft, breathable material that remains cool and dry. The brace cover is also removable and machine washable.

SECURE, NONSLIP FIT

Easily adjust the three strong fastener straps to provide a comfortable, yet secure fit. The hand and wrist immobilizer brace includes an adjustable thumb loop to separate the thumb in a neutral position.



(239) 220-5772



sales@vivehealth.com

vive
health

SR23-1R1



VIVE SPORTS COMPRESSION SOCKS

SUP3012BLK
HCPCS: A6530

Sports Compression Socks are meticulously designed to enhance blood circulation and offer superior support to your feet and ankles during physical activities. Engineered with advanced compression technology, these socks stimulate blood flow, delivering a rejuvenating experience to athletes, fitness enthusiasts, and anyone seeking comfort and performance. Crafted from premium materials, our Sports Compression Socks ensure optimal breathability and moisture-wicking properties, keeping your feet dry and odor-free. Whether you're running, cycling, or engaging in intense workouts, these socks will provide the perfect balance of compression and comfort, empowering you to push your limits and excel in your athletic endeavors.

SPECIFICATIONS

Materials:

- Spandex.
- Nylon.
- Latex-free.

Care Instructions:

- Machine Wash Cold. Hang Dry.
- Do not Bleach. Do not Iron.

Sizes:

- Small.
- Medium.
- Large.

What's included:

- One pair of compression socks

Colors:

- Black & Pantone Cool Gray C

Call now!

or send us an e-mail
now to find out more!



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com



ERGONOMIC DESIGN

Provides support to the ankle and calf muscles to minimize pressure and stiffness.



BREATHABLE MATERIAL

Moisture-wicking fabric assists with keeping feet dry and comfortable all day long.



STAYS IN PLACE

The upper-calf cuff and arch support prevent slipping during your favorite activities.

REDUCE FATIGUE AND SWELLING TO PROMOTE MUSCLE RECOVERY

The Vive Sports Compression Socks are ideal for those who play sports or exercise. With comfortable compression, these socks boost circulation and prevent muscle fatigue, for enhanced recovery.

ERGONOMIC DESIGN FOR IMPROVED CALF AND ANKLE SUPPORT

These socks offer light support for your calves and ankles, to minimize pressure and prevent stiffness in muscles and joints.

BREATHABLE FABRIC TO KEEP FEET DRY AND COMFORTABLE

Each sock features a moisture-wicking design that keeps feet dry and comfortable for all-day wear.

STAYS SECURELY IN PLACE DURING ACTIVITY

With a secure upper-calf cuff, the Vive Sports Compression Socks stay securely in place, to prevent slipping and rolling during your favorite activities.

REINFORCED ARCH SUPPORT FOR SHOCK-ABSORPTION WITH EACH STEP

Experience greater stability for the arches of your feet with reinforced support, making each step more comfortable. Great for those with plantar fasciitis.



(239) 220-5772



sales@vivehealth.com

vive
health

SR23 - 10R1

VIVE

RESTING HAND SPLINT

SUP3043

HCPCS: L3807/L3809

The Vive Resting Hand Splint is an exceptional solution designed to aid in the recovery from various hand conditions and injuries, including sprains, strains, carpal tunnel, arthritis, and trigger finger. Featuring integrated finger separators, it ensures each digit is held securely in place, promoting optimal healing and alignment. The splint combines stability and comfort, with a soft blend of foam and polyester for comfort and a removable, customizable aluminum splint for unparalleled support. Additionally, it offers adjustable straps for a personalized fit and universal sizing that fits either hand, making it a versatile choice for anyone in need of hand rehabilitation. The lightweight design and machine-washable fabric cover add to its convenience, making it an ideal companion for both healing and daily support.

SPECIFICATIONS

What's Included:

➤ 1x Resting Hand Splint

Color(s): black

Size(s): one size

Product Weight: 0.27 pounds

Product Material(s): Polyester, Foam, Aluminum

Latex free: No

Care

➤ Remove Aluminum splint

➤ Machine wash cover cold with like colors

➤ Hang dry



Call now!

or send us an e-mail
now to find out more!



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com



FINGER SEPARATORS TO PROMOTE RECOVERY

Whether you are healing from injuries like sprains or strains, or dealing with carpal tunnel, arthritis, or trigger finger, the Vive Resting Hand Splint provides dependable stability for healing and proper alignment. With integrated finger separators, this hand splint keeps each digit still and secure.



EFFECTIVE STABILITY AND SUPPORT

Ensure proper healing with support for your palm, wrist, and fingers with this brace. This stable brace prevents spasms and reduces pain, so you can focus on healing and feeling your best. The soft blend of foam and polyester creates a comfortable platform for your hand, while the removable aluminum splint ensures stability and support. The aluminum can be bent to better fit the natural shape of your hand and wrist.



ADJUSTABLE STRAP FOR PERSONALIZED COMFORT

With velcro straps, you can comfortably and securely keep your hand in a stable and steady position. The generous sizing ensures that the brace is suitable for a variety of needs, from post-operative care to joint sprains and strains.

UNIVERSAL SIZING FOR EITHER HAND

With universal sizing, this hand brace is designed to fit the left or right hand of men and women, for greater versatility and use in a variety of conditions. The lightweight design makes it perfect for transport, while the fabric cover can be machine-washed for durable and long-term use.



(239) 220-5772



sales@vivehealth.com

vive
health

SR25-2R1

VIVE STABILIZING KNEE SLEEVE

SUP3095

HCPCS: A4467

This compression knee sleeve is designed to deliver reliable support and comfort for active lifestyles or recovery needs. Featuring targeted compression, it helps improve circulation while reducing swelling, pain, and inflammation to promote faster healing after injuries or surgery. Dual spring stabilizers along the sides work together with a gel pad over the kneecap to provide cushioning and maintain proper alignment, reducing strain during movement. To ensure a secure fit, non-slip silicone grips at the top and bottom keep the sleeve in place, even during intense activities. Made from lightweight, 4-way stretch fabric, it offers breathable comfort that prevents excess heat and moisture buildup, making it suitable for all-day wear.

SPECIFICATIONS

What's Included:

➤ 1x Vive Stabilizing Knee Sleeve.

Color(s): Black

Size(s):

	Top width	Bottom width	Length
S	5.71 in	4.72 in	12.60 in
M	6.30 in	5.32 in	12.60 in
L	6.89 in	5.91 in	12.60 in
XL	7.48 in	6.50 in	12.60 in
XXL	8.07 in	6.89 in	12.60 in
XXXL	8.66 in	7.48 in	12.60 in
XXXXL	9.25 in	8.07 in	12.60 in

Material(s): 30% Nylon, 60% Rubber, 10% Polyester.

Call now!

or send us an e-mail
now to find out more!



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com



EFFECTIVE KNEE COMPRESSION

This compression knee sleeve delivers targeted pressure to improve circulation and reduce pain, swelling, and stiffness. Ideal for injury recovery, post-surgery healing, or daily joint support during activity.

STABILIZERS FOR EXTRA SUPPORT

Dual side spring stabilizers keep the knee joint steady and aligned through every movement. This added reinforcement helps prevent strain and supports active lifestyles from workouts to walking.

PROTECTIVE PATELLA CUSHIONS

The integrated gel pad surrounds the kneecap to absorb impact and provide stability. Designed to relieve pressure, it helps protect against misalignment and overuse injuries.

NON-SLIP SECURE FIT

Anti-slip silicone grips at the top and bottom prevent the sleeve from sliding down. Enjoy uninterrupted support during sports, exercise, or everyday wear without constant adjustments.

BREATHABLE ALL-DAY WEAR

Made with a 4-way stretch, moisture-wicking fabric, this knee brace offers flexible comfort. The lightweight, breathable material helps control heat and sweat for lasting wear throughout the day.



(239) 220-5772



sales@vivehealth.com

vive
health

SR25-10R1

KNEE SUPPORT

SUP3103



To see all the FAQ's,
visit [vhealth.link/jvh](https://vivehealth.link/jvh)

vive
vivehealth.com

QUICK START GUIDE

APPLICATION INSTRUCTIONS

1. Position the knee support so it is centered on your knee cap, with the support pad below it.



2. Wrap the middle strap around the back of your knee and secure it in the front.



3. Wrap the top and bottom straps around your knee, adjusting the tension for a snug, comfortable fit.



Ensure the support pad is positioned just below the kneecap for added stability, allowing your knee to move naturally without discomfort.

PROPER FIT

- Should feel snug but not restrictive.
- Adjust straps for comfort.
- Remove if numbness or tingling occurs.
- Ensure the wrap is centered over your kneecap.

QUICK START GUIDE

vive
vivehealth.com

Distributed by

vive
health
8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com

BREAK-IN PERIOD

Days 1-2	<ul style="list-style-type: none">Start with 1-2 hours, increase to 4 hours if comfortable
Days 3-7	<ul style="list-style-type: none">Gradually increase to full day as tolerated
<ul style="list-style-type: none">Take breaks every 2-3 hours to check skinRemove if pain, numbness, or irritation occursConsult doctor for nighttime use	

CARE

- Machine wash in cool water. Air dry
- Do not machine dry as this may damage the material.
- Store in cool, dry place

WARNINGS ⚠

- Consult your healthcare provider before using this product.
- Stop use and seek medical attention if pain increases or persists.
- Stop use if numbness, tingling, or discomfort occurs.

- Do not wear while sleeping unless directed by doctor.
- Do not apply too tightly - check frequently for proper fit.
- Remove immediately if swelling, skin irritation, or discomfort develops.
- This product is not intended to prevent, cure, or treat any disease.
- Not a substitute for medical treatment.
- Replace if product becomes worn or damaged.
- Check skin condition regularly during use.
- If you have sensitive skin or allergies, consult doctor before use.
- Do not wear over open wounds.
- If you have a pre-existing medical condition, consult your doctor before use.

WARRANTY

6 months warranty against defects. Contact customer service with proof of purchase.